

Spectrophotometric Determination of an Equilibrium Constant

INSTRUCTOR RESOURCES

The CCLI Initiative

Computers in Chemistry Laboratory Instruction

Learning Objective

- to determine the equilibrium constant governing the formation of $\text{Fe}(\text{SCN})^{2+}$ from iron(III) and thiocyanic acid by using the *MicroLAB Colorimeter* to measure the concentration of $\text{Fe}(\text{SCN})^{2+}$.

Procedure Overview

- colorimetric measurements for the absorbance of four different solutions of a complex ion, $[\text{Fe}(\text{SCN})^{2+}]$, are made using the *MicroLAB* Interface.
- the equilibrium constant for the formation of the complex is determined.

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT**Report Sheet****Data Table**

Complete the table while performing each trial.

	Col. A	Col. B	Col. C	Col. D	Col. E	Col. F	Col. G
Trial	Volume	Volume	Volume	Conc.	Conc.	Abs.	Conc.
1	(H ⁺)	(HSCN)	(Fe ³⁺)	[HSCN]	[Fe ³⁺]		[Fe(SCN) ²⁺]
2							
3							
4							

Stock [H⁺] (c1) =

Stock [HSCN] c2 =

Stock [Fe³⁺] c3 =

Temperature =

Provide a sample calculation used to determine the initial concentrations of [HSCN] and [Fe³⁺].

Calculations

1. Complete the "Table of Equations" with the missing equations.

TABLE OF EQUATIONS

CALCULATIONS	H ⁺	HSCN	Fe ³⁺	Fe(SCN) ²⁺
Initial Concentration	c1 ⁺	c2	c3	0
Initial Volume	V1	V2	V3	0
Initial Moles		c2V2		0
Equilibrium Moles		c2V2 - cV		cV
Equilibrium Volume	V	V	V	V
Equilibrium Concentration	(c1V1/V+c	(c2V2 / V) - c		c

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Report Sheet (page 2)

2. Provide the equations used to calculate the absorbance (from I and I₀), equilibrium concentrations and equilibrium constant.

Abs =

[Fe(SCN)₂]_{eq} =

[H⁺]_{eq} =

[HSCN]_{eq} =

[Fe³⁺]_{eq} =

K =

3. List your K values and calculate the average K and the standard deviation for your four trials.

K1	K2	K3	K4	Kave	Std. Dev.
_____	_____	_____	_____	_____	_____

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Questions/Problems

1. What effect does a dirty vial (caused by fingerprints, water spots, or lint) have on the absorbance reading for a $\text{Fe}(\text{SCN})^{2+}$ solution? Does this error cause the reported K to be too high or too low? Explain.
2. Compare the K values obtained for the different trials. Assuming a constant temperature, how would you expect the K values to compare?
3. Compare the absorbance values of your four samples. Do the values of absorbance suggest the presence of a limiting reagent in the initial solution? Explain your answer.
4. Discuss your absorbance values in relation to the concentrations of $\text{Fe}(\text{SCN})^{2+}$. Can you find a relationship? Explain.

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Questions/Problems (page 2)

5. A mixture of 2.0 ml of 2.0 *M* HCl, 5.0 ml of 0.0033 *M* HSCN and 3.0 ml of 0.0033 *M* Fe³⁺ gave a transmission current of 55 mA (*I*_t). A blank solution gave a transmission current of 82 mA (*I*₀). The calibration equation for the colorimeter used in the experiment was

$$(\text{Absorbance}) = 1400 [\text{Fe}(\text{SCN})^{2+}] + 0.011$$

Using this equation and the Table of Equations, calculate the equilibrium concentrations of Fe(SCN)²⁺, H⁺, HSCN and Fe³⁺ and then use these concentrations to calculate *K*.

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Tips and Traps

Calibration of the *MicroLAB* Colorimeter

1. Students should be instructed that the solutions of the complex are stable for a limited time; depending on the concentration. The half life, $t_{1/2}$, is approximately 1.5 hours, so they should proceed as quickly as prudent. The calibration procedure is described in the **BACKGROUND**. Known concentrations of the complex are obtained by reacting a large excess of $[\text{Fe}^{3+}]$ with HSCN. The following amounts are recommended.

2.0 M HCl (V2) ml Column A	0.0033 M HSCN(V2) ml Column B	0.20 M Fe^{3+} (V3) ml Column C	Water ml	Total Volume ml
0 (Blank)	0	2	8	10
2	1	2	5	10
2	2	2	4	10
2	3	2	3	10
2	4	2	2	10

2. The *MicroLAB* Hand Entered Spreadsheet is used to set up the calculations as follows:

Columns:	Item	Column	Item
A	Vol H^+ (v1)	B	Vol HSCN (v2)
C	Vol Fe^{3+} (v3)	D	Calc [HSCN]
E	Calc $[\text{Fe}^{3+}]$	F	Absorbance (from colorimeter)
G	Equil $[\text{Fe}(\text{SCN})^{2+}]$	H	Calc Keq.
	Average Keq.		
	Std. Dev.		

3. plot a linear regression graph.

Column F: (Abs) y-axis

Column G: $[\text{Fe}(\text{SCN})^{2+}]$ x-axis

4. The calibration equation is quite stable and can be averaged for all the boxes used.
5. The calibration equation provided in the **BACKGROUND** section works well.

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Sample Data (Page 2)

Calculations

- Complete the "Table of Equations" with the missing equations.

TABLE OF EQUATIONS

CALCULATIONS	H ⁺	HSCN	Fe ³⁺	Fe(SCN) ²⁺
Initial Concentration	c1 ⁺	c2	c3	0
Initial Volume	V1	V2	V3	0
Initial Moles	c1V1	c2V2	c3V3	0
Equilibrium Moles	c1V1+ eV	c2V2 - cV	c3V3-eV	cV
Equilibrium Volume	V	V	V	V
Equilibrium Concentration	(c1V1/V) + c	(c2V2 / V) - c	(c3V3/V) -c	c

- Provide equations used to calculate the absorbance, equilibrium concentrations and equilibrium constant in the Spreadsheet.

$$\text{Abs} = \log(I_0/I) = F$$

$$[\text{Fe}(\text{SCN})^{2+}]_{\text{eq}} = (\text{Abs} - \text{intercept}) / 2.54/\epsilon \text{ (Slope of calibration curve)} = G$$

$$[\text{H}^+]_{\text{eq}} = 4/10 + G = 0.4$$

$$[\text{HSCN}]_{\text{eq}} = D - G$$

$$[\text{Fe}^{3+}]_{\text{eq}} = E - G$$


$$K = G * (0.4 + G) / (E-G) / (D-G)$$

- List your K values and calculate the average K for your four trials.


K1	K2	K3	K4	Kave	Std. Dev.
39.69	39.70	39.63	40.41	39.86	0.37

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Sample Data (Page 2)



Sample Absorbance versus Concentration graph for FeSCN²⁺ solution



SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Suggested Answers to Questions/Problems

1. What effect does a dirty cuvette (caused by fingerprints, water spots, or lint) have on the absorbance reading for a $\text{Fe}(\text{SCN})^{2+}$ solution? Does this error cause the reported K to be too high or too low? Explain.

Less light is transmitted to the detector; this increases the absorbance reading.

The reported K will be too high. A larger absorbance indicates a higher $[\text{Fe}(\text{SCN})^{2+}]$, and also a larger K .

2. Compare the K values obtained for different trials. Assuming a constant temperature, how would you expect the K values to compare?

All the K values should be similar if the temperature does not change.

3. Compare the absorbance values of your four samples. Do the values of absorbance suggest the presence of a limiting reagent in the initial solution? Explain your answer.

Yes, the absorbance will increase if the limiting reagent increases because more $[\text{Fe}(\text{SCN})^{2+}]$ is produced (either HSCN or Fe^{3+} may become the limiting reagent).

4. Discuss your absorbance values in relation to the concentrations of $\text{Fe}(\text{SCN})^{2+}$. Can you find a relationship? Explain.

The absorbance increases if $[\text{Fe}(\text{SCN})^{2+}]$ increases, according to Beer's Law.

5. A mixture of 2.0 ml of 2.0 M HCl , 5.0 ml of 0.0033 M HSCN and 3.0 ml of 0.0033 M Fe^{3+} gave a transmission current of 55 mA (I_t). A blank solution gave a transmission current of 82 mA (I_o). If the calibration equation for the colorimeter used in the experiment was

$$(\text{Absorbance}) = 1400 [\text{Fe}(\text{SCN})^{2+}] + 0.011,$$

Use this equation and the Table of Equations to calculate the equilibrium concentrations of $\text{Fe}(\text{SCN})^{2+}$, H^+ , HSCN and Fe^{3+} and then use these concentrations to calculate K .

$$\text{Abs} = 0.17$$

$$[\text{HSCN}] = 1.53 \times 10^{-3} M$$

$$[\text{Fe}(\text{SCN})^{2+}] = 1.16 \times 10^{-4} M$$

$$[\text{Fe}^{3+}] = 8.7 \times 10^{-4} M$$

$$[\text{H}^+] = 0.40 M$$

$$K = 35 \text{ (no units)}$$

SPECTROPHOTOMETRIC DETERMINATION OF AN EQUILIBRIUM CONSTANT

Laboratory Preparation (per student station)

Equipment

- 250 ml beaker
- 10 ml Grad cylinder with stopper
- 4 plastic Beryl pipets
- cuvette
- wash bottle
- colorimeter box
- three, 25 or 50 ml burets
-

Supplies

- towel

Chemicals

Exact quantities needed are listed below. A minimum 50% excess is recommended.

- 20 ml of 0.0033 *M* KSCN
- 20 ml of 0.0033 *M* Fe(NO₃)₃ in 0.010 *M* HCl
- 10 ml of 2.0 *M* HCl
- 10 ml of 0.2 *M* Fe(NO₃)₃ in 0.010 *M* HCl for calibration

Safety and Disposal

- no special precautions necessary