

## Recrystallization

7/22/2022

1

### Recrystallization technique

---

**Recrystallization** is a laboratory technique used to purify solids based on their **different solubilities**.

A **small amount** of solvent is added to a flask containing an impure solid. The contents of the flask are **heated** until the solid dissolves.

Next, the solution is **cooled**.

A more pure solid precipitates, leaving impurities dissolved in the solvent.

Vacuum **filtration** is used to isolate the crystals.

The waste solution is discarded. The impurity either remains dissolved or removed before cooling by absorption on **decolorizing carbon** ( Norit)

7/22/2022

2

## Selection of solvent: characteristics

---

An ideal crystallization solvent should have:

- a) It should dissolve a large amount of the substance at its boiling point and only a small amount at room temperature.
- b) It should dissolve the impurity at low temperature or not at all.
- c) Its boiling point should be lower than the melting point of the solid substance.
- d) Upon cooling, the solvent should yield well formed crystals.
- e) It should not react with the solute.
- f) It should be safe (nonflammable, nontoxic)
- g) It should be inexpensive

## Selecting a solvent:

---

Carry out solubility tests on the compound to be re-crystallized.

Place 0.1 g of the powdered solid into a clean test tube. Add the solvent dropwise with shaking.

If the **solid dissolves**, the solvent is **not fit** for crystallization.

If it does **not** dissolve, **warm** the solution gently to boiling. If it partially dissolves, add more solvent and warm.

**Cool** the solution and observe the **appearance** of crystals.

**Note:** if the solid does not dissolve in a large amount of the hot solvent, then it is not fit for re-crystallization.

## Summary of Recrystallization Steps

---

The **objective** is to dissolve the solute in a **minimum** amount of the hot solvent.

1. Select the **proper** solvent ( see the notes below)
2. Place a **small** amount of solid in an Erlenmeyer flask. Add a boiling chip.
3. Add a **small quantity** of the appropriate hot solvent .
4. **Heat** to boiling on a steam bath (hot plate or Bunsen burner can be used).
5. To the boiling solution, add more **hot solvent** in small portions with stirring.
6. Continue addition of solvent until all solute dissolves at boiling.
7. If colored impurities are present, add a small amount of **charcoal** to the solution (not boiling solution) and boil.
8. Filter the hot solution by gravity using a preheated short stem funnel. Apply heat to dissolve the solid.
9. Cool the solution to crystallize the product. Slower cooling may lead to a higher purity product, so it's common practice to allow the solution **to cool to room temperature** before setting the flask in an ice bath.

7/22/2022

5

## Summary of Recrystallization Steps-Continued:

---

10. Crystals usually begin forming on the bottom of the flask.

\* It is possible to aid crystallization by **scratching** the flask with a glass rod at **the air-solvent junction** (assuming you are willing to purposely scratch your glassware).

The scratch increases the glass surface area, providing a roughened surface on which the solid can crystallize.

- Another technique is to '**seed**' the solution by adding a small crystal of the desired pure solid to the cooled solution. Be sure the solution is cool, or else the crystal could dissolve.
- If no crystals fall out of solution, it's possible too much solvent was used. Allow some of the **solvent to evaporate**. If crystals do not spontaneously form, reheat/cool the solution.

11. Use vacuum filtration to collect the pure crystals

12- Wash the crystals with a small amount of the cold solvent and press out the solvent.

13. The crystals can be air dried or placed in a vacuum desiccators.

7/22/2022

6

## Decolorizing:

If there are colored impurities, the solution (in step 3 above) is cooled.


Charcoal ( 2-3% of the weight of sample is added to the cold solution which is then boiled with stirring and filtered. ( Steps 4 & 5 above).

7/22/2022

7

## Summary of steps

1)Transfer the solid 2) Heat Solvent 3)Pour hot solvent gradually 4)Pure solid precipitates  
Impurities remain in solution 5)Filter




7/22/2022

8

## Preparation and Re-crystallization of Acetanilide:

---

Acetanilide is a white solid that can be prepared by reacting aniline and acetic anhydride:



**Aniline + Acetic Anhydride      Acetanilide + Acetic Acid**

Basic Information:

Compound	Molecular Mass	M.P. /(B.P.)	Density
Acetanilide	135.17	114 °C	-----
Aniline	93.13	(184 °C)	1.022 g/mL
Acetic Anhydride	102.09	(138 °C)	1.082 g/mL

7/22/2022

9

## Experimental Procedure:

---

- 1- Weigh 3 grams of crude acetanilide in a 250 mL-Erlenmeyer flask.
- 2- Add boiling water gradually to dissolve all the white solid.
- 3- Cool a bit and add **0.2 g** of decolorizing carbon ( Charcoal).
- 4- Boil the solution again, filter out charcoal using a preheated Buchner funnel and suction flask. Then while hot, transfer the filtrate to clean beaker.
- 5- Cool the filtrate to room temperature then in ice (crystals should appear).
- 6- Use suction filtration to collect the crystals on a Buchner funnel. Wash the crystals with a small amount of cold water.
- 7- Keep the suction on to squeeze out all the solvent. Air dry the crystals.
- 8- Weigh the dry product, calculate the percentage yield and determine its melting point.
- 9- Collect the product in a paper and write your name and submit it to your instructor.

7/22/2022

10