

Experiment 1

Melting Points

7/19/2022

1

Melting Points:

Physical properties are often utilized by chemists in **identifying** an organic compound. These physical properties include :
color, odor, physical state, melting point (M.P.), boiling point (B.P.), density (d), infrared (IR) spectrum, nuclear magnetic (NMR) spectrum and ultraviolet (UV) spectrum.

As long as the physical constants are determined under **standard conditions** (temperature, pressure, etc.), they are invariant and, therefore, useful in helping to determine the **identity of unknown substances**.

There are a number of **reference books** that contain tables of physical properties and physical constants of compounds.

One of the most common is the *Handbook of Chemistry and Physics*. If the physical properties of an unknown compound are identical to the physical properties of a compound listed in the tables, the two compounds are **probably** the same.

7/19/2022

2

Definitions:

The **melting point of a solid** is defined as the **temperature** at which the liquid and solid phases are in equilibrium.

The **freezing point of a liquid** is the same temperature as the melting point of its solid. However, freezing points are **rarely** measured in practice because they are more difficult to determine. One reason for this is that solidification may not occur at the correct temperature due to the phenomenon of **supercooling**.

Supercooling occurs when a liquid is cooled below its freezing point does not solidify.

Thus, in practice, most **melting points** are determined as capillary melting points, which can be done quickly with a small amount of sample in a capillary tube. A **capillary melting point** is defined as the temperature range over which a small amount of solid in a thin walled capillary tube first visibly softens (first drop of liquid) and then completely liquefies.

7/19/2022

3

Melting Points and Mixed Melting points:

Melting points are determined for three reasons.

- 1-If the compound is a **known one**: the melting point will help to **characterize** the sample in hand.
- 2-If the compound is **new** then the melting point is recorded in order to allow **future characterization** by others.
- 3- The **range** of the melting point is indicative of the **purity of the compound**; an impure compound will melt over a wide range of temperatures.

*Pure organic compounds generally have **sharp melting points**. An **impurity** lowers the melting point and widens the range. **Impurities** in a solid cause a melting point **depression** because the impurity **disrupts the crystal lattice energies***

A technique for proving the identity of an unknown compound is the **mixed melting point**. Advantage is taken of the depression of melting points of mixtures to prove whether two compounds having the same melting points are identical.

If X and Y are identical, then a mixture of the two will have the same melting point; but if X and Y are not identical, then a small amount of X in Y or of Y in X will cause the melting point to be lowered.

7/19/2022

4

Eutectic Mixtures:

Eutectic mixture is defined as a **mixture** of two or more components which usually **do not interact** to form a new chemical compound but, which at certain ratios, inhibit the crystallization process of one another resulting in a system having a **lower melting point** than either of the components

7/19/2022

5

Eutectic Mixtures:

The melting point behavior of impure compounds is best understood by consideration of a simple binary mixture of compounds X and Y (Fig. 1).

This melting point-composition diagram shows the melting point behavior as a function of composition. The **melting point of a pure** compound is the temperature at which the **vapor pressures** of the solid and liquid are equal.

But in dealing with a mixture the situation is different. Consider the case of a mixture of 75% X and 25% Y.

At a temperature below ET, the eutectic temperature, the mixture is solid Y and solid X.


At the eutectic temperature the solid begins to melt. The melt is a solution of Y dissolved in liquid X. The vapor pressure of the solution of X and Y together is less than that of pure X at the melting point; therefore, the temperature at which X will melt is lower when mixed with Y.

As the temperature is raised, more and more of solid X melts until it is all gone at point M (temperature m). The melting point range is thus from ET to m.

7/19/2022

6

Eutectic Mixtures



7/19/2022

7

The Thiele apparatus

This apparatus achieves stirring and **uniform heat distribution** by convection. It is filled to the base of the neck with silicone oil (the oil expands on heating). and equipped with a thermometer clamp to hold the thermometer.


The tube is heated at the **base of the bend**.

The bulb of the thermometer should be **halfway** down the tube to assure uniform heating.

Melting points are also easily determined in a beaker.

The beaker can be heated on a hot plate or a Bunsen burner.

Do not discard the oil used in the apparatus because it will be necessary to determine a number of melting points in future experiments.



7/19/2022

8

Filling Melting Point Capillaries

1-The dry sample is ground to a **fine powder** on a watch glass or a piece of paper on a hard surface using the flat portion of a spatula.

2-It is formed into a small pile and the melting point **capillary forced down** into the **pile**.

3-The sample is **shaken into** the closed end of the capillary by rapping sharply on a hard surface or by dropping it down a 50 cm length of glass tubing onto a hard surface.

Note: The height of the sample should be no more **than 2-3 mm**.

7/19/2022

9

Determining the Melting Point

The accuracy of the melting point depends on the accuracy of the thermometer, so the first exercise in this experiment will be to calibrate the thermometer.

1- **Melting points of pure**, known compounds will be determined and deviations recorded so that a correction can be applied to future melting points.

2- The most critical factor in determining an accurate melting point is the **rate of heating**. At the melting point the temperature rise should not be greater than 1°C per minute. This may seem extraordinarily slow, but it is necessary in order that heat from the bath be transferred **equally** to the **sample** and to the **glass and mercury of the thermometer**.

7/19/2022

10

Part 1: Calibration of the Thermometer

Determine the melting point of standard substances (Table 1) over the temperature range of interest.

Include both temperatures making up the **melting point range**. Use your data to graph a correction curve for your thermometer.

Compound	Observed m.p. Range (Ti-Tf)	Average	Literature m.p.
1,4-Dichlorobenzene			
Acetanilide			
Salicylic Acid			
4-Nitrobenzoic acid			

7/19/2022


11

Part 2: Melting Points of Unknown

To determine the melting point of an unknown provided by your instructor:

Prepare **two capillaries of each unknown**. Run a very **fast** determination on the first sample to ascertain the approximate melting point and then cool the melting point bath to just below the melting point and make a **slow**, careful determination using the other capillary.

You may wish to run a **mixed melting** point as well.



7/19/2022

12